

SKYLA VB1

CHEMISTRY ANALYSER



The Skyla VB1 Clinical Chemistry Analyser is a compact system designed to perform blood biochemistry tests – fast, reliably and with high accuracy.

Its easy-to-use and simple operation requires no specialised training and enables you to obtain results in about 15 minutes with the precision and accuracy equivalent to those of a clinical laboratory analyser.

Features:

- ▶ Fast and reliable clinical chemistry analysis
- ▶ Compact, bench top unit
- ▶ 200µl sample volume
- ▶ Can use whole blood, serum or plasma
- ▶ Species specific reference ranges
- ▶ Stores up to 50,000 test results
- ▶ Built-in printer
- ▶ 2 year warranty
- ▶ Large LCD touch
- ▶ Almost no maintenance required



1. Insert whole blood into test disk



2. Place test disk on drawer



3. Get results

