## THE QUICKVET ANALYSER

COAGULATION, FIBRINOGEN & BLOOD TYPING



## Bringing lab-on-a-chip technology to your practice.

The QuickVet Analyser provides analytical results on a par with large laboratory analysers for PT, aPTT, fibrinogen and blood typing.

The state of the art 'Lab-on-a-Chip' analysis system profits from minimised manual intervention, increased speed of analysis, increased precision of data, significant cost reduction, less sample consumption and minimised operator exposure to hazardous material.

## The QuickVet Analyser is:

- Easy to use
- ▶ Accurate, reliable and cost effective
- Uses only 30ul of blood sample
- ➤ Tests results are easy to understand and are available in 2-4 minutes for coagulation, 14 mins for fibrinogen and 5 mins for blood typing.
- Available tests include:
- ► Coagulation PT/aPTT Combo Test
- Equine Fibrinogen The first in-clinic and stableside Fibrinogen test available
- Canine Blood Typing Canine DEA 1.1 Test
- Canine Fibrinogen
- ▶ Feline Blood Typing





