



vepalabs
VETERINARY PATHOLOGY

INSIGHT HCT METER



The Insight HCT meter is a low cost handheld meter for the measurement of canine and feline haematocrit and total haemoglobin in capillary and venous whole blood.

The Insight HCT meter utilises optical reflectance for determination of haematocrit (HCT) and total haemoglobin. The appropriate species chip (cat or dog) and test strip are inserted into the meter and a drop of whole blood is applied to the test strip. The meter displays haematocrit (HCT) and total haemoglobin results in just 30 seconds, reporting HCT as a percentage (%).



Features

- ▶ Accurate and precise
- ▶ Handheld 30 second test time
- ▶ One drop whole blood (Venous, Capillary or EDTA)
- ▶ Stores 250 results
- ▶ HCT results reported as %
- ▶ Room temperature strip storage
- ▶ LCD Screen