# i-SMART 30 VET



**ELECTROLYTE ANALYSER** 

Fast, accurate and maintenance free in-practice electrolyte analysis

#### **Fast & Accurate**

- ▶ Gold Standard ISE electrolyte analysis
- ► Tests electrolytes Na\*, K\* and CI\*
- ▶ 60ul sample volume
- Results in 35 seconds
- ▶ Specific animal species reference ranges

## **All-in-One Cartridge**

- ▶ No user maintenance
- ► Auto-calibration
- Cartridge contains electrodes and all consumables

## **Features & Advantages**

- ▶ 300,000 test data storage
- ► Large touch screen LCD
- Veterinary specific software
- ▶ Battery or mains powered



#### **Measured Electrolytes**

Measuring Range
20 ~ 250mmol/L
05 ~ 20.0 mmol/L
20 ~ 250 mmol/L

