

# G-PET PLUS+

## VETERINARY GLUCOMETER



*The new and improved  
veterinary glucometer*

## g-PET Plus+ makes testing quick, easy and economical

The g-PET Plus+ glucose veterinary monitoring system features:

- ▶ Biosensor auto code technology
- ▶ 1.5 microlitre blood sample
- ▶ 4 reminder alarms for worry free testing
- ▶ 8 second countdown
- ▶ 300 result memory with date and time
- ▶ Accurate results
- ▶ Big LCD screen
- ▶ Hypoglycaemic and hyperglycaemic indicator alarms
- ▶ PC download option



## Testing made easy with just 3 step

**Step 1:** Insert the test strip (metre will turn on automatically)

**Step 2:** Apply blood

**Step 3:** Results in 10 seconds



**STEP 1**



**STEP 2**



**STEP 3**

