

URIT CA 11

HEMATOLOGY CALIBRATORS

CALIBRATOR

Assay values and ranges provided by URIT Medical
Model skyla SH50



2024/11/25

Instrument	Parameter	LOT	24071404M3
SH50	WBC (WOC) × 10 ⁹ /L	9.92	±0.40
	WBC (WIC) × 10 ⁹ /L	9.62	±0.40
	RBC × 10 ¹² /L	3.84	±0.20
	HGB g/dL	8.4	±0.3
	HGB g/L	84	±3
	MCV fL	62.5	±3.0
	PLT × 10 ⁹ /L	219	±20

Note:

1. Take these values for reference only when the calibrator within the validity period. Do not use it as the calibrator expired.
2. Do not take the above values a reference as the lot number or reagent system being changed.
3. Shake the calibrator up and down for 1-2 minutes as taking it out from refrigerator and place it in room temperature (15°C~30°C) for 15 minutes, and then doing the test.

URIT Medical Electronic Co., Ltd.

No. D-07 Information Industry District, High-Tech Zone, Guilin, Guangxi 541004, P.R.China.

Tel: +86(773)2288586 Fax: +86(773)2288560

Web: www.urit.com Email: service@uritest.com

EC-REPRESENTATIVE:

Shanghai International Holding Corp. GmbH (Europe)

Eiffestrasse 80, 20537 Hamburg, Germany

Manufactured for Shenzhen URIT under contract